

PS-2076

Product Description:


P-DERM® PS-2076 is a 0.22mm thickness silicone gel and acrylic trilaminate consisting of silicone gel skin contact adhesive, polyurethane and medical pressure sensitive acrylic. PS-2076 is supplied with a 75 micron polycarbonate liner on the silicone gel adhesive side and an 84# bleached kraft siliconized liner on the acrylic side of the construction.

Construction / Properties:

Property	Value	Test Method
Total Thickness	0.22mm	QSP-726
Silicone Gel Coat Weight	150 gsm	QSP-724
Polyurethane Carrier	25 micron	QSP-724
Acrylic Adhesive Coat Weight	37 gsm	QSP-724
Silicone Gel Adhesion	0.5 N/25mm	QSP-723
Acrylic Adhesion	4.3 N/25mm	QSP-722
Continuous Use Conditions	10 - 70 C	QSP-754
Biocompatibility	Cytotoxicity: Pass Primary Dermal: Pass Buehler Sensitization: Pass	ISO 10993
Recommended Sterilization Method	Ethylene Oxide	--

Features:

- Excellent Instantaneous Tack
- Atraumatic Removal from Skin
- Hypoallergenic
- Ability to be Removed and Repositioned
- Excellent Converting Properties



Specific tests should be performed by the end user to determine the product stability for the particular application.

For Additional Information:

E-mail: sales@polymerscience.com
Toll Free: +1 888.533.7004
Web: www.polymerscience.com
Revision: 010319